

## Problems

### 1. Images do not open when clicking the Image icon.




#### a. Solution

- Exit SkyVue.
- Close report window.
- Reopen report and click the image icon .


### 2. Images not loading quickly

- a. Load times for historical images can take up to 5 minutes to load.
- b. SkyVue will open and the user will see a green status bar in the lower right indicating that SkyVue is transferring and loading the images.
- c. A team has been engaged in the migration of historical images to a new server and anticipates completion by late summer 2022.




### 3. CINE-when done viewing, the user clicks X to close view. The screen goes blank.

#### o Solution

- **CINE**-when done viewing, **click the CINE icon**  again.
- **Do NOT click X**, which closes the viewer and results in a blank screen.

### 4. Customizing Shortcut tools

- a. **Right click image.** Depending upon the image, shortcut tools will vary.
- b. Click the **wrench icon** to customize tools.
- c. Select the desired tool from Available Tools.
- d. Click the arrow to move to the right Displayed Tools area.
  - o Click OK.



### 5. Printing Images from SkyVue

## SkyVue Troubleshooting for Providers, Nurses, Clinical Staff

### Cerner PowerChart, Firstnet and Ambulatory EDUCATION

- The print option does not work.
- Solution
  - Use the Snipping Tool to capture an image, copy it to a Word document, and print.

#### 6. To view other films on the same patient:

- Solution: Add and default the History Toolbar.
  - Click **Workspace** at the top of the page.
  - Select **Workspace Gadgets**.
  - Click **History**.
  - Images are selectable from the history toolbar at the bottom of the page.
  - Left click the **History** tab, drag and drop next to the Thumbnail tab.
  - Click **Workspace**.
  - **Click Save Workspace.**

