

Pyxis MedStation: How to do a Global Facility Patient Search

1. Log in to the Pyxis MedStation
2. Select All Available Patients

The screenshot shows the Pyxis MedStation interface. At the top, the user is logged in as Beeman, Sheryl. The main area displays a list of patients under the heading 'All Available Patients'. The list includes columns for Name, Visit ID, and Location. A red box highlights the 'All Available Patients' tab and the search bar. A text box explains that if a patient is not found in the 'All Available Patients' list, the 'Facility Search' button should be used.

Name	Visit ID	Location
ATLASQES, MUNSON A 46 y (8/12/1970) M	1710300017	
ATLASTEST, GOOFY 61 y (10/15/1955) M	1719900001	
ATLASTEST, GOOFY 61 y (10/15/1955) M	1719900003	OPM
CERNBABYTEST, TOPAZ 5 m (2/13/2017) F	1704400009	NICU 1100/A
CERNNAPLUS, SEVEN 50 y (9/22/1966) M	1702400005	B4 4019/A
CERNPORTALTEST, APPLE 47 y (1/1/1970) F	1720800003	EDM
CERNPORTALTEST, ORANGE 47 y (6/7/1970) F	1720900001	EDM
CERNPOUTEST, ACDC 73 y (5/17/1944) M	1721300016	GFS GFS/11
CERNPOUTEST, AEROSMITH 54 y (11/8/1962) M	1721300018	GFS GFS/13
CERNPOUTEST, BACH 45 y (8/13/1971) M	1721400007	IR IR/10
CERNPOUTEST, BEYONCE 29 y (4/18/1988) F	1721400009	PG-Mat 1113/A
CERNPOUTEST, CCR 51 y (3/19/1966) M	1721300020	B2S B2S/10
CERNPOUTEST, CHER 33 y (5/16/1984) F	1721400011	PG-Mat 1115/A

3. Enter Last Name or visit ID and the search button

The screenshot shows the Pyxis MedStation interface. The search bar is highlighted with a red box, and the text 'zzzzcer' is entered. The search button is also highlighted with a red box. Below the search bar, the text 'There are no results to display.' is shown.

Note: You cannot add these types of patients to your "My Patients" list until they are transferred in the system to your nursing unit.