


Respiratory Specimen Collection Device

Unless otherwise stated in notes, multiple tests can be collected in ONE device

WARNING: Specimens collected in 1 ml of media, with more than one test requested, risk **NOT** having enough for COVID testing

WARNING: Specimens collected in 3 ml of media cannot have a rapid flu test performed



Device	Picture of media	Used for flu test?	Used for Respiratory Panel?	Used for CoVID testing?	Notes:	Media storage pre collection/ media storage post collection
NP swab collected and placed in paper wrapper (done with point of care devices)	NA	Yes	No	No		Room Temp/Room Temp
NP swab collected in Unitranz media		Yes	Yes	Yes	If 3 test are required, 2 Unitranz with NP collection required	Room Temp/ Refrigerated Transport to lab as quickly as possible

Respiratory Specimen Collection Device

Unless otherwise stated in notes, multiple tests can be collected in ONE device

WARNING: Specimens collected in 1 ml of media, with more than one test requested, risk **NOT** having enough for COVID testing

WARNING: Specimens collected in 3 ml of media cannot have a rapid flu test performed



Device	Picture of media	Used for flu test?	Used for Respiratory Panel?	Used for CoVID testing?	Notes:	Media storage pre collection/ media storage post collection
NP swab collected in CoPan UTM 1 ml collection device		Yes	Yes	Yes	If 3 test are required, 2 CoPan UTM with NP collection required	Room Temp/ Refrigerated Transport to lab as quickly as possible
NP swab collected in CoPan UTM 3 ml collection device		No	Yes	Yes		Room Temp/ Refrigerated Transport to lab as quickly as possible

Respiratory Specimen Collection Device

Unless otherwise stated in notes, multiple tests can be collected in ONE device

WARNING: Specimens collected in 1 ml of media, with more than one test requested, risk **NOT** having enough for COVID testing

WARNING: Specimens collected in 3 ml of media cannot have a rapid flu test performed

Device	Picture of media	Used for flu test?	Used for Respiratory Panel?	Used for CoVID testing?	Notes:	Media storage pre collection/ media storage post collection
NP Swab collected in BD Universal Transfer media		No	Yes	Yes		Room Temp/ Refrigerated Transport to lab as quickly as possible
NP Swab in M5 Viral Transport media		No	No	Yes		Refrigerated/ Refrigerated Transport to lab as quickly as possible