

Depression Screening Recommendation Workflow		Screening Diagnosis Code	
1	Obtain PHQ-9 paper form from patient (if used)	ICD-10	Description
2	Go to the Vital Signs component on the Ambulatory Nursing Workflow	Z13.31	Depression Screening
3	Click the down arrow next to Vital Signs and select the Intake PowerForm	CPT-4	Description
4	Navigate to the PHQ2 and PHQ-9 section of the Intake PowerForm	96127	Depression Screening
5	Enter responses for PHQ2 and/or full PHQ-9 (based on practice requirement) Note: If the PHQ-2 Score is 1 or more, the PHQ-9 must be completed	G0444	Depression Screen (Medicare Only)
6	If the patient answers anything other than Not at all to Thoughts Better Off Dead or Hurting Self (Question #9 on the PHQ-9 Screening), notify the provider	HCPCS	Description
7	Time spent administering the depression screening tool must be documented in the PHQ-9 Comments field	G8510	Depression Screen Negative no plan reqd
8	Continue documenting the patient intake	G8431	Depression Screen Positive w/ Plan
9	On the Problems and Diagnosis section of the Intake form, if the patient does not have an active chronic diagnosis of depression, bipolar disorder and/or dysthymia to select this visit, add Z13.31 Depression Screening to Diagnosis (Problem) being Addressed this Visit	G8511	Depression Screen Positive w/o plan
10	When the Intake is completed, click the green checkmark to sign		
11	Go to Recommendations component on Ambulatory Workflow		
12	Locate Depression Screening recommendation and click on Orders		
13	Select appropriate Depression Screening CPT-4 code: Brief Emotional/Behavioral Assessment or Depression Screen (Medicare Only)		
14	Ordering Physician box: select Order and communication type Cosign Required and enter the ordering provider's name		
15	Click on Depression Screening to open the details. Click Actions and select the appropriate HCPCS code: Depression Screen negative or Depression Screen Positive w/ or w/o Plan		
16	Click Orders for Signature button and associate the Z13.31 Depression Screening diagnosis and click Sign		